# ELSODENT

2 x 5 ml

Adhésif universel  
Universal bonding

## Healbond Max

G-PHARMA - 35, avenue du Gros Chêne - B.P. 10279 HERBLAY  
95617 CERGY-PONTOISE Cedex FRANCE  
elsodent@elsodent.com - www.elsodent.com

### Restaurations indirectes

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GB

### Description

Two-component universal adhesive system including a PRIMER and an ADHESIVE. The Healbond Max Primer provides effective etching to enamel and dentin without the need for an etching gel, thus simplifying the bonding procedure. The material's chemistry allows for compatibility with all self-cure or dual-cure resin cements and core build-up materials. Dentists can therefore use Healbond Max for direct and indirect procedures, without additional bonding system. HEMA-free formulation ensures a durable bonding.

### Indications for use

**Note** : Healbond Max is compatible with all composites, materials use to cast inlay cores, self-cure, dual-cure and light-cure resin-cements and core buildup materials, and it does not require the use of any dual cure Activator.

### Direct Bonding Application

- Light-cured composite restorations.
- Composite and fractured ceramic/metal crowns repairs.
- Sealing of hypersensitive and/or exposed root surfaces.
- Core build-up materials (self-cured, light-cured, or dual-cured).

**Indirect Bonding Applications** (with self-cure, dual-cure and light-cure resin-cements and light-cure composites for veneers only)

- Veneers.
- Indirect Inlay cores, and in conjunction with self-cure, dual-cure, and light-cure resin cements).
- Zirconia, ceramic (lithium disilicate-based and alumina-based), composite, and metal-based (precious and non-precious) inlays, onlays, crowns, bridges.
- Endodontic posts.

### Safety notes

#### Warnings

- Healbond Max is a corrosive and flammable liquid which contains ethanol and uncured methacrylate. These components may cause pulp damage, skin, eyes and mucosa irritation. They can lead to contact dermatitis with persons at risk.
- Avoid eyes contact, to prevent any irritation or damage of the cornea. In case of eyes contact, rinse immediately and thoroughly with copious amount of water. Seek for a medical advice.
- Avoid skin contact to prevent any irritation or potential allergic reaction. In case of skin contact, remove the product immediately with a cotton pellet and thoroughly wash and rinse with soap and water. If any allergic reaction appeared (erythema), seek for medical advice.
- Avoid any contact with mucosa to prevent inflammatory reactions. In case of contact, remove the product immediately. Once the restoration is done, rinse thoroughly and ask the patient not to swallow the water. If any allergic reaction of the mucosa persists, seek for medical advice.

#### Precautions

- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.

### Contraindications

Healbond Max is contraindicated in patients with a history of allergic reactions to methacrylate resins or any other of the components. Do not use Healbond Max for direct pulp capping in case of pulp exposure.

### Clinical procedures

Use one disposable applicator brush for Healbond Max Primer and another one for Healbond Max Adhesive.

### Direct Restorations

#### 1. Bonding Composites to dentin/enamel

- Prepare the cavity, using a rubber dam to isolate the oral cavity. Clean unprepared tooth structure with a fluoride-free cleaning paste; rinse thoroughly with water spray and air dry.
- Apply Healbond Max Primer to the enamel/dentin surface using a disposable applicator brush. Scrub the dentin and enamel surface with a brushing motion for 20 seconds.
- Air dry until no fluid movement is visible.
- Shake Healbond Max Adhesive bottle briefly. Apply to the enamel/dentin surface with light brushing motion for 15 seconds.
- Gently air dry to get a uniform bond film.
- Light-cure for 10 seconds.
- Apply composite per manufacturer's instructions for use.

#### 2. Post Bonding and Core Buildup

- Prepare the surface of the post per manufacturer's instructions.
- Apply Healbond Max Primer in the root canal preparation and on all the cavity walls involved in the bonding procedure, then Air dry until no fluid movement is visible.
- Apply Healbond Max Adhesive in the root canal preparation and all the walls of the cavity clobbered by the bonding procedure.
- Gently air dry to get a uniform bond film. Excesses of Healbond Max Adhesive in the canal can be removed with a dry applicator brush or an absorbent paper point.
- Light-cure for 10 seconds.
- Apply Healbond Max Primer in the post surface and air dry until no fluid movement is visible.
- Apply a thin coat of Healbond Max Adhesive on the post surface. Gently air dry to get a uniform bond film.
- Light-cure for 10 seconds.
- Apply the resin cement of your choice onto the post, seat the post, vibrate the post slightly avoiding the possibility of trapped air.
- Lighth-cure the resin-cement per the manufacturer's instructions. Use the core material of your choice. All self-cure and dual-cure core materials can be built in bulk as light-cured Healbond Max is compatible with ail self-cure and dual-cure core buildup material.

### Important notices

The use of Healbond Max Adhesive is not necessary if INNOCEM is the resin-cement of your choice. Healbond Max gets high enamel bond strengths without the use of etching. However, etching gel can be used on enamel and/or dentin if a total etch or selective enamel etch technique is desired. In that case, after step 1, etch Enamel or Dentin surfaces for 15 seconds, then rinse thoroughly. Proceed to step 2.

### Indirect restorations

- Internal surface preparation (Ceramic, metal, Zircon, alumine, composite)**
  - Abrase the surface with 50µm alumina with a pressure of about 30 psi (0.2 MPa).
  - For ceramic restorations only: etch the restoration with hydrofluoric acid for 1 minute. Rinse thoroughly and well dry with air.
  - Apply Healbond Max Primer on the internal surface of the restoration and air dry until no fluid movement is visible.
  - Apply a coat of Healbond Max Adhesive on the internal surface of the restoration. Gently air dry to get a uniform bond film.
- Light-cure for 10 seconds.

### Important notices

The use of Healbond Max Adhesive is not necessary if INNOCEM is the resin-cement of your choice. Healbond Max Primer can be used as silane replacement.

### 2. Tooth Surface Preparation

- Clean the preparation, Wash thoroughly with water spray and air dry.
- Apply Healbond Max Primer to the enamel/dentin surface. Scrub the surface with a brushing motion for 20 seconds.
- Air dry until no fluid movement is visible.
- Shake Healbond Max Adhesive bottle briefly. Apply the material to the enamel/dentin surface with light brushing motion for 15 seconds.
- Gently air dry to get a uniform bond film.
- Light cure for 10 seconds.

### Important notices

The use of Healbond Max Adhesive is not necessary if INNOCEM is the resin-cement of your choice. Healbond Max produces high enamel bond strengths without the use of etching. However, etching gel can be used on enamel and/or dentin if a total etch or selective enamel etch technique is desired. In that case, after step 1, etch Enamel or Dentin surfaces for 15 seconds, then rinse thoroughly ensuring that all etchant has been removed. Proceed to step 2.

### 3. Cement Application and Seating of Restoration

#### Note

Healbond Max is compatible with all type of inlay cores, composites, self-cure, dual-cure and light-cure resin-cements and core buildup materials. It does not require the use of any dual cure Activator.

#### a. Veneers

Choose the proper shade of cement. Apply the resin cement of your choice to the veneer, seat the veneer, and cure the cement per manufacturer's instructions for use.

#### b. Crowns, bridges, inlays, onlays and metal-based restorations

Apply the resin cement of your choice to the restoration or preparation, seat the restoration, and cure the cement per manufacturer's instructions for use.

### Storage & shelf-life

- Healbond Max should be tightly closed immediately after use.
- Do not expose the product to sun-light or warm sources. Store at a temperature between 2°C and 8°C.
- Do not use after the expiration date printed on the box.

### Presentation & references

BMAX-10 – Universal bonding, 2 x 5ml bottles (Primer + Adhesive).

### Important notice

- Do not throw away the box and the instruction for use until the product is finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.

### Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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