# HARMONY

GB

## Description

Harmony vinyl polysiloxane impression materials have been designed to make high precision dental impressions. The extreme hydrocompatibility of the material, its mat appearance and its thixotropy, makes it one of the most performing impression materials available, in term of clinical features.

## Times

- **HLF-100** : Low viscosity, fast setting (orange).
  - Working time : 1 minute 30 seconds
  - Total setting time at 23°C : 2 minutes

- **HLN-100** : Low viscosity, normal setting (green).
  - Working time : 2 minutes 15 seconds
  - Total setting time : 3 minutes 30 seconds

- **HM-100** : Medium viscosity, (purple).
  - Working time : 1 minute 45 seconds
  - Total setting time : 3 minutes 10 seconds

- **HPN-600** : Putty soft, normal setting (yellow).
  - Working time : 2 minutes 30 seconds
  - Total setting time : 3 minutes 30 seconds.

- **HPF-600** : Putty soft, fast setting (blue).
  - Working time : 1 minute 45 seconds
  - Total setting time : 2 minutes

## Contraindication

Harmony is contraindicated in patients with a history of allergic reaction to vinylpolysiloxane materials, or any other of the components.

## Instructions for use

- Remove the cap from the cartridge and check that the two exits are not obstructed.

- Attach a mixing tip (yellow ring type) and an intraoral tip (yellow), on the cartridge.

- Place the cartridge in the mixing gun.

- When using the cartridge for the first time, bleed a small amount of material onto a mixing pad to ensure that both base and catalyst will be extruded equally. Discard the purged amount.

## Impression procedure

### One-step technique

*Recommended materials: HPN-600 and HM-100*

- Gently remove the retraction cord, or wash off retraction material, from the sulcus of the prepared teeth. Flush and dry.

- Prepare the putty material by manually mixing the base and catalyst pastes (50/50) and load a perforated metallic impression tray.

- While the tray is being loaded, apply the wash material around the cleaned and dried preparation.

- Depending on the number of teeth involved, initiate application during the putty preparation, so the tray loading and application of material around the preparation are completed at the same time.

- Be careful not to exceed the working time of the material to avoid incomplete or inaccurate impressions.

- Slowly place the loaded tray in the mouth, and hold it in position passively.

- Allow the impression material to fully set, then remove the impression from the mouth.

- Wash the impression thoroughly and place it in disinfectant solution, prior to send it to the technician laboratory.

- The duration should comply with the manufacturer IFU of the solution.

### Two-step technique ("wash technique")

*Recommended materials: HPF-600 and HLF-100 or HLN-100*

- Prepare the putty material by manually mixing the base and catalyst pastes (50/50) and load a perforated metallic impression tray.

- Gently remove the retraction cord, or wash off retraction material, from the sulcus of the prepared teeth. Flush and dry.

- Slowly place the loaded tray in the mouth, and hold it in position.

- Allow the impression material to fully set, then remove the impression from the mouth.

- Excise undercuts and silicone present in the inter-dental space, and cut spillways, if required.

- Clean and dry the impression thoroughly to provide for optimal adhesion between the putty and the light materials.

- Insert the light impression material into the preliminary impression at the appropriate sites and apply around the cleaned and dried preparations, while keeping the intra-oral tip immersed in the paste, to avoid the formation of bubbles.

- Replace the preliminary impression into the mouth and hold it in position passively.

- Allow the impression material to fully set, then remove the impression from the mouth.

- Wash the impression thoroughly and place it in disinfectant solution, prior to send it to the technician laboratory.

- The duration should comply with the manufacturer IFU of the solution.

## Storage & shelf-life

Store in a dry, cool and well ventilated place, away from humidity. Store at ambient temperature not exceeding 25°C. When stored under such conditions, Harmony has a shelf-life of 24 months.

## Presentations & references

- **Harmony light & medium** : 2 x 50ml dual cartridge + 10 yellow mixing tips

- **HM-50** (medium consistency – normal setting)

- **HLF-50** (light consistency – fast setting)

- **HLN-50** (light consistency – normal setting)

- **Harmony putty soft** : 2 x 300ml jars (base + catalyst)

- **HPF-600** (putty soft – fast setting)

- **HPN-600** (putty soft – normal setting)

## Important notice

Keep the cartridge with the mixing tip in place. Never re-use the cap of the cartridge! Because re-using the original cap to close the cartridge could induce a cross-contamination and lead to a premature polymerization of the product. Do not touch the impression material with latex gloves. They can inhibit its polymerization. Do not throw away the box and the instruction for use until the product is completely finished. High humidity could damage the packaging on which important instructions and legal information are printed. Storage at a too high temperature could cause to age prematurely and inhibit its polymerization. Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

## Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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