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Self-etching, dual-cure permanent resin cement

1 x 5ml dual syringe +10 mixing tips + 10 intra-oral tips

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ELSODENT

Description/Indication

Dual -curing, self-etching, self-adhesive cement, universal use, for definitive cementation of all types of indirect restorations, whatever their composition. The extra-fine intraoral tips are an answer to all clinical situations.

InnoCem can be used either as a "Total Etch" or a "Self-Etch" procedure. The additional etching of the tooth structure with G Etch increases the adhesive strength.

Indications

Definitive cementation of ceramic, Zirconia, composite or metal crowns and bridges, on natural teeth or implants. In combination with Healbond MP, universal bonding, InnoCem is indicated for :

- Definitive cementation of metal, ceramic, or composite inlays, onlays and veneers.
- Definitive cementation of Maryland bridges.
- Definitive cementation of endodontic metal or fiber posts.

Safety notes**Warnings**

- InnoCem contains polymerizable methacrylates which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosapersists, seek medical attention.

Precautions

- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.

Contraindication

InnoCem is contraindicated in patients with a history of allergic reaction to methacrylate resins or any other of the components.

Instructions for use**Tooth structures preparation**

- Clean the tooth structures, rinse with a water spray, and lightly air dry. Do not over dry. The tooth structures should be just dry enough that the surface has a slightly glossy appearance.
- For veneers and Maryland bridges, tooth structure must be etched.

Treatment of metal, composite, Zirconia and aluminum oxide ceramic restorations

- Blast or have the composite, metal, Zirconium or non etchable ceramic restoration surfaces to be luted, blasted with aluminum oxide. Clean the blasted surfaces with alcohol and dry with air free of water and oil.
- If the restorations are onlays, inlays or veneers, apply Healbond MP and wait 15 seconds before drying slightly for 5/7seconds until the solvent has evaporated. Do not light cure

Treatment of etchable ceramic restorations

- Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.
- Rinse thoroughly with water for 15 seconds and dry with air free of water and oil.
- Apply Healbond MP following IFU and wait 15 seconds before drying slightly for 5/7seconds until the solvent has evaporated. Do not light cure.

Preparation of Root Canals&Treatment of glass fiber-reinforced posts

- Choose and try the selected post, in the canal.
- Clean the root canal with a solution of sodium hypochlorite (NaOCl). Rinse immediately and thoroughly with water, dry with paper tips.
- Clean the post with alcohol and dry it with air free of water and oil. When using Elsdent Fiber Post, application of a silane is not required. Other fiber-reinforced posts should be treated in accordance with their instructions for use.
- Apply Healbond MP on both the post and canal walls, following IFU and wait 15 seconds before drying slightly for 5/7seconds until the solvent has evaporated. Do not light cure.

Cementation

To avoid pulpal irritation on vital teeth, cover deep areas, close to the pulp with an anti-bacterial liner (Healdent), or a calcium hydroxide liner.

1. When using a new syringe of InnoCem, remove dual syringe cap, dispense and discard a small amount of material. Be sure the material is flowing freely from both ports.
2. Attach a mixing tip on the dual syringe.

Cementation of indirect restorations (inlays, onlays, crowns and bridges)

1. Remove the provisional restoration.
2. Remove temporary cement residus of the prepared tooth or implant abutment surface.
3. Rinse thoroughly.
4. Dry without excess the prepared tooth or implant abutment.
5. Try in restoration to verify its good adaptation, especially at the limits of the preparation.
6. Depress syringe plungers to begin the flow of material. Do not force. Extrude enough product to deposit a thin uniform layer in the inside of the definitive restoration.
7. Seat the restoration carefully onto the prepared tooth or implant abutment, leaving overflow cement.
8. Light-cure the excesses of cement on the edges, for 3 to 5 seconds until it will turn to a gel. Take them off with a probe.
9. For composite and ceramic restorations, light-cure the cement through the restoration.
10. Hold the restoration firmly in place for 2minutes 30seconds max. (self-curing mode).
11. Polish the joint with strips or rubber discs cups.
12. Check occlusion.

Note

In case of low abutment devitalized tooth, an etching will increase the adhesive strength.

Cementation of a post in a root canal**Option 1**

Use Healbond MP in conjunction with Core D, core build up material (see IFU of these products).

Option 2

1. Attach an Xfine intra-oral tip to the mixing tip for direct application in the root canal.
2. Insert the tip as deeply as possible in the root canal and apply InnoCem, beginning apically. Keep the intraoral tip in the cement and slowly move upwards.
3. Place the post in the root canal filled with cement and apply moderate pressure to hold the post in position and rotate the post slightly during insertion to avoid the inclusion of air bubbles.
4. Let the cement flowing out of the canal.
5. Light cure InnoCem for 5 seconds to remove the excesses.

Cementation of veneers and Maryland bridges

1. Apply Healbond MP on the etched surfaces of the tooth, lightly rubbing the dentin surface for 15 seconds Dry slightly for 5 seconds until the solvent has evaporated, and light-cure for 10 seconds.
2. Inject directly InnoCem into the veneer or Maryland bridge surface. Gently place the element on the tooth leaving the cement overflow.
3. Cure through the veneer in place, on the vestibular face, and distance to margins, using the "boost power" of the Led curing light, for about 5 seconds. For Maryland bridge, cure at the limits of the element for 5 seconds.
4. After having removed cement excess, light-cure the entire surface, as well as the joint, for 10 seconds.

Setting time

Total setting time is 3 minutes and 10 seconds, in the self-curing mode at 37°C, and 20 seconds (translucent shade) to 30 seconds (opaque shades) in the light-curing mode.

Storage and shelf-life

Do not expose the product to sun-light or warm sources. Store at a temperature not exceeding 25°C. When stored under such conditions, InnoCem has a shelf-life of 24 months.

Presentation & references

IC-10-TA2 : 1 x 5ml syringe + 10 mixing tips + 10 intra-oral tips – Translucent-A2 shade

IC-10-OA2 : 1 x 5ml syringe + 10 mixing tips + 10 intra-oral tips – Opaque A2 shade

IC-10-OA3 : 1 x 5ml syringe + 10 mixing tips + 10 intra-oral tips – Opaque A3 shade

Important notice

Do not throw away the box and the instruction for use until the product is completely finished. High humidity could damage the packaging on which important instructions and legal information are printed. Storage at a too high temperature could cause to age prematurely and inhibit its polymerization. Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the instructions for use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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