Notre Dam

GB

Definition

Notre Dam is a methacrylate-based light-cure barrier (liquid dam), used in dental procedures. It is hard when set, however easy to remove from the gum, dentin walls, undercuts and embrasures.

Bleaching

- To protect gum during bleaching sessions at the chair.

Endodontic

- To fill and hermetically seal the gap between the dam and the tooth during endodontic sessions.
- To built a temporary missing wall of the tooth, to keep the irrigation liquid in the operative field.

Prosthesis

- To block undercuts, in case of receding gum, before impression.

Safety notes

Warnings

- Notre Dam contains polymerizable methacrylate resins which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

Precautions

- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.

Contraindications

Notre Dam is contraindicated in patients with a history of allergic reaction to methacrylate resins or any other of the components. Do not use Notre Dam in cases of pulpular exposure, pulpitis.

Instruction for use (bleaching)

- Remove the cap from the luer lock syringe and secure a disposable applicator.
- Rinse thoroughly the designed area and slightly dry the surface of the teeth, gingival third, the embrasures or the gaps to fill.
- Apply Notre Dam onto the gingival and the embrasures. If some product flows in excess onto the enamel, cure the excess for 5 seconds and use a mouth spatula to remove it. For the embrasures, be sure that the resin injected from the vestibular side will appear on the lingual side. This extension will add retention to the resin barrier.
- When the barrier is applied for bleaching, extend the resin barrier beyond the last tooth to be bleached.
- The resin barrier should be about 5-6mm height and 2 mm thick.
- Light-cure for 20 seconds.

Clean-up procedure

- Rinse thoroughly the treated area to eliminate the whitening gel and vacuum.
- Lift the polymerized dam, using a probe or a spatula. The interproximal dam is eliminated by using a probe and/or floss.
- Rinse again the gum and interproximal areas.

Storage & shelf-life

Do not expose the product to sun-light or warm sources. Store at a temperature not exceeding 25°C. When stored under such conditions, Notre Dam has a shelf-life of 36 months.

Presentation & reference

Blue colour for a good contrast.
ND - 4 - 2 x 1,2ml syringes + 10 disposable applicators

Important notice

- Do not throw away the box and the instruction for use until the product is completely finished.
- High humidity could damage the packaging on which

- important instructions and legal information are printed.
- Storage at a too high temperature could cause a degree prematurely and inhibit its polymerization
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and its use for any purpose not explicitly stated in the instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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