Definition

Light-cure opaquer used to mask :

- Exposed metal surfaces during the repair of defective ceramic or composite veneers.
- Dark dentin during direct class III or IV restoration procedures.

Safety notes**Warnings**

- Opaq contains polymerizable methacrylate resins which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

Precautions

- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.

Contraindications

Opaq is contraindicated in patients with a history of allergic reaction to methacrylate resins or any other of the components. Do not use Opaq in cases of pulpular exposition, pulpitis.

Instructions for use**Restoration of fractured porcelain teeth**

- Provide adequate isolation to the operative field by using a rubber dam.
- Roughen the metal surface by using an intraoral blasting device or a diamond bur.
- Bevel the defective area of the ceramic or composite material with a diamond bur.
- Dry with oil-free air.
- Etch, the ceramic or composite surface and the exposed metal surface, for 10-15 seconds (for ceramic, use 9.6% hydro-fluoric acid gel, following manufacturer's recommendations). Rinse thoroughly with water and dry with oil-free air.
- Apply a silane agent on the ceramic or composite surface. Wait 60 seconds and dry with oil-free air. A contact between the silane agent and the metal surface does not negatively influence the bonding performance.
- Apply a layer of 0.5mm thickness of Opaq, on the conditioned metal surface exclusively. Polymerize for 20 seconds. If necessary apply a second opaquer layer and polymerize it.
- Apply a bonding agent on the entire surface to be repaired (opaquer and ceramic/composite material), blow it to a thin layer, and polymerize for 20 seconds.
- For proximal restorations, use a transparent matrix band and wedges.
- Fill the cavity by applying a composite in the chosen VITA shades, in layer of 2mm. Polymerize each layer for 20 seconds at a light intensity of at least 500 mW/cm² or 10 seconds at a light intensity of at least 1100mW/cm².
- Finish and polish the restoration.

To mask dark dentin under composites restorations

- Prepare the cavity in conventional manners.
- Etch the dentin (15-20 seconds) and enamel (40 seconds). Rinse and slightly dry.
- Apply a layer of 0.5mm thickness of Opaq. Polymerize for 20 seconds.
- Apply a bonding agent on enamel, Opaq and free-dentin. Slightly dry and photo-polymerize for 20 seconds.
- For proximal restorations, use a transparent matrix band and wedges.
- Fill the cavity by applying a composite in the chosen VITA shades, in layer of 2mm. Polymerize each layer for 20 seconds.
- Finish and polish the restoration.

Storage & shelf-life

Do not expose the product to sun-light or warm sources. Store in a cool, dry place, at a temperature not exceeding 25°C when stored under such conditions, Opaq has a shelf-life of 36 months.

Presentations & reference

Available in A3 VITA shade.
OPA-2 - 1 x 1,2ml syringe + 10 disposable applicators

Important notice

- Do not throw away the box and the instruction for use until the product is completely finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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