Description

Resin-based temporary luting cement. Well-balanced adhesive properties provide good prosthesis retention. Versatile in use : as a temporary cement on natural abutment teeth or, as a definitive cement for crowns on implants. Its translucent beige shade makes Seal Temp the perfect material for temporary cementation of definitive aesthetic elements (veneers, onlays). Excellent resistance to oral environment and mechanical strength assure good long-term performance. Low film thickness of cement.

Safety notes**Warnings**

- Seal Temp contains polymerizable methacrylate resins which may irritate to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

Precautions

- Use of rubber dam or adequate isolation is recommended.
- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.

Contraindication

Seal Temp is contraindicated in patients with a history of allergic reaction to methacrylate resins or any other of the components.

Instructions for use**Seal Temp, double syringe form**

- Remove dual syringe cap, dispense and discard a small amount of material. Be sure the material is flowing freely from both ports.
- Install a mixing tip on the dual syringe
- Gently depress syringe plungers, to begin the flow of material. Do not use excessive force. Extrude enough material to make a thin layer of cement inside the temporary crown.
- Hold the temporary crown firmly in place.
- Trim excesses immediately, with a probe or a spatula, before polymerization occurs.
- Keep the dual syringe with the mixing tip in place until the next use.

Seal Temp, 2 syringes form

- Dispense approximately equal amount of Part A and Part B pastes onto a mixing pad.
- Mix the pastes thoroughly with a spatula for 10 to 15 seconds
- Spread a thin layer of cement inside the temporary crown and hold the crown firmly in place.
- Trim excesses immediately.

Notice

Seal-Temp is retentive enough for low height abutment teeth.

Storage & shelf-life

Do not expose the product to sun-light or warm sources. Store at a temperature not exceeding 25°C. When stored under such conditions, Seal Temp has a shelf-life of 36 months.

Presentation & references

- STD-32** - 2 x 10ml double syringes + 10 mixing tips + 10 intra-oral tips. (to be used with a dispensing gun)
STD-16 - 2 x 5ml double syringes + 20 mixing tips.
STD-8 - 1 x 5ml double syringe + 10 mixing tips.
ST-25 - 2 x 12,5g syringes (base + catalyst).

Important notice

Do not throw away the box and the instruction for use until the product is completely finished. High humidity could damage the packaging on which important instructions and legal information are printed. Storage at a too high temperature could cause to age prematurely and inhibit its polymerization. Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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