CORE D & CORE D Flow

GB

Definition

Dual-cure, highly resistant composites for core build-up. Core D & Core D Flow are available in a tooth-colored shade which is ideal in situations where esthetics of the core is important. Permanent cementation of endodontic fiber-posts, in conjunction with Healbond MP, universal self-etching adhesive system.

Physical data/Characteristics of the material

- Working time : 90 secondes
- Intra oral setting time + working time (self-curing mode) : 3 minutes 10 secondes
- Flexural strength : 120 Mpa
- Compressive strength : 300 Mpa
- Radioopaque

Safety notes

Warnings

- Contains polymerizable methacrylate resins which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

Precautions

- The use of this product is restricted to dental professionals.
- The product must not be used for any purpose other than specified in the indications.
- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.

Contraindication

The product is contraindicated in patients with a history of allergic reaction to methacrylate resins.

Interactions

Cement containing Eugenol (ZnO cements) may inhibit the polymerization of Core D and Core D Flow and leads to discoloration.

Composition

Hydrophilic di-methacrylate resins - modified urethane di-methacrylate resins - barium aluminosilicate glass - amorphous silica - Camphorquinone (CQ) - accelerators - stabilizers - pigments.

Clinical procedures

First use of a new cartridge

- Remove and discard the cap of the cartridge. Under no circumstances should the cap of the cartridge be used to reseal the cartridge as this could result in the inclusion of air bubbles, inducing a cross-contamination (base/catalyst) and leading to a premature polymerization of the product.
- Check whether both cartridge openings are free and remove obstructing material, if any, with a suitable instrument.
- Before installing a mixing tip, bleed the cartridge onto a mixing pad until the base and catalyst pastes are at the same level.

Important notice

When Healbond MP is used with Core D or Core D Flow, there is no need to use a dual cure activator.

Core built up

- Prepare the cavity in the usual manners.
- Refer to manufacturer's instructions for tooth conditioning and adhesive application. For optimum

bonding, we recommend the use of Healbond MP, self-etching universal adhesive system. Apply Healbond MP without activator and light-cure for 10 seconds.

- Once the surfaces have been properly treated, they must be kept uncontaminated.
- Install a mixing tip + an intra-oral tip on the cartridge (Core D), or on the dual-barreled syringe (Core D, Core D Flow).
- Gently depress the trigger of the dispensing gun, or the syringe plungers, to begin the flow of material. Do not use excessive force.
- Apply Core D or Core D Flow, directly in the treated cavity. Core D can be applied without a matrix in place because it does not flow. A banded matrix retainer or clear crown form is needed for Core D Flow.
- Light-cure for 30 seconds at a light intensity of at least 600 mW/cm² or 15 seconds at a light intensity of at least 1100mW/cm².
- Once has set, use rotary instruments to complete the crown preparation (final shaping and finishing) in the same clinical session.
- After use, leave the mixing tip in place: it will protect the nozzle until the next application.

Cementation of endodontic fiber-post

- Refer to manufacturer's instructions for canal walls conditioning and adhesive application. For optimum bonding, we recommend the use of Healbond MP (self-etching universal adhesive system) without activator. If core build-up will immediately follow post cementation, treat all available tooth structure beyond post space with Healbond MP for core build-up bonding.
- Treat surface of post according to the post manufacturer's instructions.
- Once the surfaces have been properly treated, they must be kept uncontaminated.
- Apply Healbond MP without activator and light cure for 10 seconds.
- Install a mixing tip + an intra-oral tip on the cartridge (Core D), or on the dual-barreled syringe (Core D, Core D Flow).
- Gently depress the trigger of the dispensing gun, or the syringe plungers, to begin the flow of material. Do not use excessive force.
- Apply the material directly onto the treated post surface and into the treated post hole preparation.
- Immediately place post into final position and stabilize.
- Once stabilized, light-cure, all accessible areas of the post for 20 seconds, at a light intensity of at least 600 mW/cm²

10. Proceed with the core build-up.

Storage & shelf-life

- The product must be used by the expiration date indicated on the package.
- Do not expose the product to sun-light or warm sources.
- Store at a temperature between 5°C and 25°C. When the product is stored in the refrigerator, it should be brought back to room temperature 15 minutes minimum before use, to restore its normal viscosity and curing time.
- Use the product at room temperature.

Presentations

Core D is available in dual cartridge and dual syringe. Core D Flow is available in dual syringe.

Important notice

- Leave the used mixing tip in place on the cartridge until next use.
- Do not throw away the box and the instruction for use until the product is finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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