

CORE D

Dual-cure core build-up material

Formulation without addition of HEMA, TEGDMA

ADVANTAGES

Very thixotropic. Does not flow. No need for a matrix. Improves operating field vision

Hardness equivalent to dentin. No different feeling while burring

Radio-opaque

High flexural strength > **140 MPa**

High compressive strength ~ **330 MPa**

High conversion rate

INDICATIONS

Core build-up on non-vital teeth

Sealing root canal posts associated with a dual-curing bonding type HEALBOND MP



Very thixotropic



REFERENCES & PRESENTATIONS

CDB-25	25ml cartridge + 20 mixing tips + 20 intra-oral tips - white shade
CDD-25	25ml cartridge + 20 mixing tips + 20 intra-oral tips - dentin shade
EJ-100	100 yellow mixing tips
IJ-100	100 yellow intra-oral tips
IV-50	50 green intra-oral tips

MD : Class IIa certified by Tüv Rheinland (0197), excl. tips, Class 1.