

# CORE D FLOW

Dual-cure core build-up fluid material  
Formulation without addition of HEMA, TEGDMA

## ADVANTAGES

Flowable, with thin tips, to allow injection in small and narrow cavities.

**Hardness equivalent to dentin.** No different feeling while burring

High conversion rate = 68% (1h)

High flexural strength = **136 MPa**

High compressive strength ~ **330 MPa**



## INDICATIONS

Core build-up on non-vital teeth

Sealing root canal posts associated with a dual-curing bonding type HEALBOND MP

## REFERENCES & PRESENTATIONS

<b>CDF-5</b>	2 x 2,5ml syringes + 10 mixing tips + 10 intra-oral tips - dentine
<b>EMP-100</b>	100 brown mixing tip
<b>IMP-100</b>	100 intra-oral tips

*MD : Class IIa certified by Tüv Rheinland (0197), excl. tips, Class 1.*

## IMPORTANT NOTICE !

ELSODENT formulations of **CORE D products family** allows **HEALBOND MP** to become **DUAL** without using of **an activator**.