

# HARMONY

Vinyl polyxilane impression materials

## ADVANTAGES

Specially developed to provide a very precise and highly hydrocompatible impression material

**Thixotropic.** Doesn't start to flow until the Putty is in place. Improves patient comfort. Accurate and "easy to take" impression.

Matte appearance for a perfect and fast reading. Correction possible right away, if needed

High compressive strength

Good dimensional stability < 0.2%

Excellent coverage after distortion > 99%



## VISCOSITY & SETTING TIME

Low viscosity, regular setting  
2-steps impression (Wash Technique)

Low viscosity, fast setting  
2-steps impression (Wash Technique)

Medium viscosity, regular setting.  
For "one step" impression procedure

Putty Soft, regular setting

Putty Soft, fast setting

## REFERENCES & PRESENTATIONS

<b>HPN-600</b>	Putty soft, regular setting 2 x 300ml (base + catalyst)
<b>HPF-600</b>	Putty soft, fast setting 2 x 300ml (base and catalyst)
<b>HM-100</b>	Medium viscosity, regular setting - 2 x 50ml cartridges + 10 mixing tips
<b>HLN-100</b>	Low viscosity, regular setting - 2 x 50ml cartridges + 10 mixing tips
<b>HLF-100</b>	Low viscosity, fast setting - 2 x 50ml cartridges + 10 mixing tips
<b>HM-200</b>	Medium viscosity, regular setting - 4 x 50ml cartridges
<b>HLN-200</b>	Low viscosity, regular setting - 4 x 50ml cartridges
<b>HLF-200</b>	Low viscosity, fast setting - 4 x 50ml cartridges
<b>EJ-100</b>	100 yellow mixing tips
<b>IJ-100</b>	100 Intra-oral tips

MD : Class I