

ELSOTEMP

Eugenol-free, fluoridated temporary filling cement

ADVANTAGES

Because of its unique plasticizer, no expansion will occur during the setting process

No pulp or gum irritations (PH=7)

Fast-setting

Waterproof

Does not dry in the pot, even when opened

Eugenol-free: does not hinder polymerization of resins and composites

Contains fluoride: cariostatic action

Soft before setting to allow a perfect application during the filling process

Hard after setting

Stick to the walls of the tooth: remains in mouth 1 or 2 months if needed

Easy to remove



CONSISTENCY

Grey	soft
White	standard
Pink	hard

REFERENCES & PRESENTATIONS

TPS-B	Temporary filling cement 3 x 28g white
TPS-R	Temporary filling cement 3 x 28g pink
TPP-B	Temporary filling cement 6 x 28g white
TPP-G	Temporary filling cement 6 x 28g grey
TPP-R	Temporary filling cement 6 x 28g pink

MD : Class IIa certified by Tüv Rheinland (0197).